EN Sphygmomanometer

ES Esfigmomanómetro

DE Sphygmomanometer

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
EN English

Introduction

Thank you for choosing Rossmax Palm Type Aneroid Sphygmomanometer. Blood pressure measurements determined with Palm Type Aneroid Sphygmomanometer are equivalent to those obtained by a trained observer using cuff/stethoscope auscultation method.

Palm Type Aneroid Sphygmomanometer is protected against manufacturing defects by an established International Warranty Program. For warranty information, you can contact the manufacturer, Rossmax International Ltd. or your local distributors.

Attention: Consult the accompanying documents. Please read this manual carefully before use. Please be sure to keep this manual.

Preliminary Remarks

Palm Type Aneroid Sphygmomanometer complies with European Regulation underlying the Medical Device Directive and bears the CE mark "CE1639". The quality of device has been verified and conforms to the provisions of the EC council directive 93/42/EEC (Medical Device Directive) Annex I essential requirements and applied harmonized standards..

EN ISO 81060-1: 2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type.

To ensure accurate measurements, this monitor is recommended to be re-calibrated every two years.

About Blood Pressure

It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk.

For reliable monitoring and reference of blood pressure, keeping long term records is recommended. Please download the blood pressure log (chart for recording and tracking crucial info: date, time, blood pressure, pulse rate and body condition) at www.rossmax.com.

Applying the Cuff

1. Place the cuff over the bare upper arm, wrap it with the tube pointing your palm, and the artery mark over your main artery.
2. The edge of the cuff should be at approximately 1.5 to 2.5 cm above the inner side of the elbow joint.
3. Using the correct cuff size is important for an accurate reading.

Measurement Procedures

1. Place the stethoscope head over the main artery, underneath the artery mark of the cuff. (Fig. ①)
2. With the valve closed, press the bulb and continue pumping to a value 20-30 mmHg above your normal blood pressure. (Fig. ②)
3. Open the valve to deflate the cuff gradually at a rate of 2-3 mmHg per second. (Fig. ③)
4. Record the onset of Korotkoff sound as the systolic pressure, and the disappearance of these sounds as diastolic pressure.

5. After measurement is completed, fully open the valve to release any remaining air in the cuff.
Note: 1. Sit upright in a chair, and take 5-6 deep breaths. Avoid leaning back while the measurement is being taken.

2. Do not cross the legs while sitting and keep the feet flat on the floor during measurement.

Cautionary Notes

1. The unit contains high-precision assemblies. Therefore, avoid extreme temperatures, humidity, and direct sunlight. Avoid dropping or strongly shocking the main unit, and protect it from dust. 2. The unit should not be operated by children to avoid hazardous situations.

3. The maintenance of the Palm Type Aneroid Sphygmomanometer must be performed by an authorized service. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact your local distributor or Rossmax International Ltd.

4. To stop operation at any time, open the air release valve, and the air in the cuff will be rapidly exhausted.

5. The lot number is for service tracking.

6. During the blood pressure measurement, blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.

7. Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.

8. Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.

Storage and Maintenance

1. After measurement, wrap cuff around sphygmomanometer and bulb and store in the carrying case together with the stethoscope.

2. Always deflate cuff completely before storage.

3. Never touch the cuff fabric or parts with a sharp instrument.

4. Sphygmomanometer accuracy can be checked visually: simply be certain the needle rests under "0" when the unit is fully deflated. Under normal condition, it is recommended that calibration be completed every two years.

5. Clean the cuff with a slightly damp sponge or cloth. If necessary, it may be washed in cool water with mild soap. Remove bladder prior to washing. Rinse and air dry. Do not iron. Clean the bladder and bulb with a slightly damp sponge or cloth. To protect from cracking or discoloring, wipe with a soft cloth moistened with ethanol. Clean valve and gauge with a soft dry cloth.

Safety Information

1. Rossmax Palm Type Aneroid Sphygmomanometer is exclusively designed for blood pressure measurement at the upper arm or thigh on healthy skin.

2. Do not inflate to more than 300 mmHg.

3. Do not take your blood pressure for more than 3 minutes.

4. Wait for at least 5 minutes before repeated measurements.

Specifications

Measurement Range: 0-300 mmHg
Accuracy: ±3 mmHg
Operation environment: 10°C~40°C; 85% RH max
Storage environment: -10°C~60°C; 10%~90% RH
Dimensions: 163(L)x60(W)x45(H)mm
Weight: 170±5g

* Specifications are subject to change without notice.

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www.rossmax.com

Warranty Card









This instrument is covered by a 2 year guarantee from the date of purchase, cuff and accessories are not included. The guarantee is valid only on presentation of the warranty card completed or stamped by the seller/dealer confirming date of purchase or the receipt. Opening or altering the instrument invalidates the guarantee. The guarantee does not cover damage, accidents or non-compliance with the instruction manual. Please contact your local seller/dealer or buying source or www.rossmax.com.

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