EN Blood Pressure Monitor

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S1 3'AAA size, 1.5V batteries included in the carrying case with the unit.

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EN English

Introduction

Blood pressure measurements determined with S150 are equivalent to those obtained by a trained observer using cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard, Electronic or Automated Sphygmomanometers. This unit is to be used by adult consumers in a home environment. The patient is an intended operator. Do not use this device on infants or neonates. S150 is protected against manufacturing defects by a one year warranty period. For warranty information, you can contact the manufacturer, Rossmann International Ltd, or your local distributor.

Attention: Consult the accompanying documents. Please read this manual carefully before use. For specific information on your own blood pressure, contact your physician.

Please be sure to keep this manual.

How This Unit Works

S150 uses the oscillometric method to detect your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure, which is equivalent to the air pressure. The measurement of your blood pressure is based on this baseline pressure. After the cuff inflates to block your blood in the artery, the deflation process starts. During the deflation of the wrist cuff, the monitor is detecting the pressure oscillations generated by the beat-to-beat pulse. After detecting the pressure oscillations during this period of time will cause measurement error. After detecting the amplitude and the slope of the pressure oscillations during the deflation process, your S150 will determine for you the systolic and diastolic pressures, and your pulse rate is detected at the same time.

Preliminary Results

The blood pressure monitor complies with the European regulations and bears the CE mark. The quality of the device has been verified and conforms to the provisions of the EC council directive 93/42/EEC (Medical Device Directive).

Annex I: General requirements and applied harmonized standards.

EN 1060-1: 1995/A2: 2009 Non-invasive sphygmomanometers - Part 1 - General requirements

EN 1060-3: 1997/A3: 2009 Non-invasive sphygmomanometers - Part 3 - Supplementary requirements for electro-mechanical blood pressure measuring systems

EN 1060-4: 2004 Non-invasive sphygmomanometers - Part 4: Test Procedures to determine that overall system accuracy of non-invasive sphygmomanometers

ISO 8160-2: 2013 Non-invasive sphygmomanometers - Part 2: Clinical investigation of automated measurement type.

This blood pressure monitor was designed for long service time. To ensure accurate measurements, this monitor is recommended to be re-calibrated every two years.

Blood Pressure Standard

The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, classifying blood pressure ranges into 4 stages. (Ref. The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure-Complete Report JNC-7, 2003). This blood pressure monitor is based on the classification of blood pressure ranges into 4 stages.

Display Explanations

EE / Measurement Error: Measure again. Wrap the cuff correctly and keep wrist steady during measurement. If the error keeps occurring, return the device to your local distributor or service center.

E1 / Critical Artery Abnormality: Measure again. If the error keeps occurring, return the device to your local distributor or service center.

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E3 / Data Error: Remove and replace the batteries. If the error keeps occurring, return the device to your local distributor or service center.

E4 / Measurement Range: Measure again. If the error keeps occurring, return the device to your local distributor or service center.

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Movement Detection

The "Movement Detection" helps reminding the user to remain still and is indicating any body movement during measurement. The specified icon appears once a body movement has been detected during and after each measurement.

Note: It is strongly recommended that you consult your physician if the IHB icon appears often.

Installing Batteries

Open the battery cover in the right direction. Install battery in the right position. Battery Type: 2 AAA size. Replace the cover and click in the other end to secure the battery cover.

Replacing the Cuff

Clamp the cuff to left wrist with palm facing up as Ø.

Make sure the edge of the cuff is about 1 cm from the palm as Ø.

In order to ensure accurate measurements, fasten the velcro strap securely around your wrist so there is no extra space between the cuff and the wrist as Ø. If the cuff is not wrapped tight enough, the measurement values may be false.

Do not service or maintain device and cuff while in use.

Applying the Cuff

1. Clean the skin of the patient's arm with a cloth. Dry the skin with a clean cloth.

2. Apply the cuff to the patient's arm.

3. Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous shunt.

4. In order to ensure accurate measurements, fasten the velcro strap securely around your wrist so there is no extra space between the cuff and the wrist as Ø. If the cuff is not wrapped tight enough, the measurement values may be false.

5. If your physician has diagnosed you with poor circulation on your left arm, carefully place the cuff around your right wrist as shown in Ø.

Correct Measuring Posture

1. Place your elbow on a table so that the cuff is at the same level as your heart as Ø. Note: Your heart is located slightly below your armpit a bit to the left of the middle of your chest. Relax your entire body, especially the area between your elbow and fingers.

2. If the cuff is not at the same level as your heart or if you can not keep your arm completely still throughout the reading, use a soft object such as a folded towel to support your arm as Ø. Do not allow hard objects to come in contact with the wrist cuff.

3. Turn your palm upwards.

4. Sit upright in a chair, and take 5-6 deep breaths. Avoid leaning back while the measurement is being taken as Ø.


Measurement Procedures

Important Notes:

• These are a few helpful tips to help you obtain more accurate readings.

• Blood pressure changes with every heartbeat and is in constant fluctuation throughout the day.

• Blood pressure recording can be affected by the position of the user, his or her physiologic condition and other factors. For greatest accuracy, wait one hour

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EN 2'AAA size, 1.5V batteries included in the carrying case with the unit.

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