EN English Introduction

Thank you for choosing Rossmax GA series, the shock and splash resistant backlighted digital sphygmomanometer. Blood pressure measurements determined with GA series are equivalent to those obtained by a trained observer using cuff/stethoscope auscultation method.

GA series is protected against manufacturing defects by an established International Warranty Program. For warranty information, you can contact the manufacturer, Rossmax International Ltd. or your local distributors.

Attention: Consult the accompanying documents. Please read this manual carefully before use. Please be sure to keep this manual.

Preliminary Remarks

GA series complies with the European regulations and bears the CE mark "CE 1639". The quality of the device has been verified and conforms to the provisions of the EC council directive 93/42/EEC (Medical Device Directive), Annex I essential requirements and applied harmonized standards.

EN 1060-1: 1995/A2: 2009 Non-invasive sphygmomanometers - Part 1 - General requirements
EN 1060-3: 1997/A2: 2009 Non-invasive sphygmomanometers - Part 3 - Supplementary requirements for electro-mechanical blood pressure measuring systems

Blood Pressure Standard

The National High Blood Pressure Education Program Coordinating Committee has developed a blood pressure standard, classifying blood pressure ranges into 4 stages.

(Ref. The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure-Complete Report JNC-7, 2003).

This blood pressure classification are based on historical data, and may not be directly applicable to any particular patient. It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk before making any healthcare decision based on the device's readings.

For reliable monitoring and reference of blood pressure, keeping long-term records is recommended. Please download the blood pressure log (chart for recording and tracking crucial info: date, time, blood pressure, pulse rate and body condition) at www.rossmax.com.

Name/Function of Each Part

- On/Off/Zeroing key
- Marker key
- Air Tube and Connector
- Inflating Bulb
- Air Release Valve
- Arm Cuff

Display Explanations

- Heartbeat Mark
- Weak Battery
- Inflation Pressure (during measurement)
 - Systolic Pressure (after measurement)
- Deflation Rate (during measurement)
 - Diastolic Pressure (after measurement)
- Error Mark
- Calibration Reminder
- Rotate the clip through 90°.
- Loosen the screw on the battery cover.
- Install the batteries with correct polarities.
- Replace the battery cover and tighten the screw.
- Rotate the clip back.
- You need to replace the batteries when:
 - Low battery icon appears on display.
 - The On/Off/Zeroing key is pressed and nothing appears on display.

Caution:
 1. Batteries are hazardous waste. Do not dispose them together with the household garbage.
 2. There are no user serviceable parts inside. Batteries or damage from old batteries are not covered by warranty.
 3. Use exclusively brand batteries. Always replace with new batteries together. Use batteries of the same brand and same type.

Applying the Cuff

- Plug the connector into the gauge.
- Place the cuff over the bare upper arm, wrap it with the tube pointing your palm, and the artery mark over your main artery.
- The edge of the cuff should be at approximately 1.5 to 2.5 cm above the inner side of the elbow joint. If the index line falls within the range of the arm circumference indicator, the cuff circumference is suitable, otherwise you may need a cuff with a different circumference.
- Using the correct cuff size is important for an accurate reading.
- To stop operation at any time, open the air release valve, and the air in the cuff will be rapidly exhausted.

Measurement Procedures

- Place the stethoscope head over the main artery, underneath the artery mark of the cuff. (Fig. ①)
- Press the On/Off/Zeroing key. After zeroing, GA series is ready to measure. (Fig. ②)
- With the valve closed, press the bulb and continue pumping to a value 30-40 mmHg above your normal blood pressure (Fig. ③). The back light will be turned on once the pressure is inflated over 30 mmHg, and throughout the measurement process.
- Open the valve to deflate the cuff gradually at a rate of 2-3 mmHg per second (The deflation rate is displayed on LCD) (Fig. ④).
- Press the Marker key to record the onset of Korotkoff sound as the systolic pressure, and press the Marker key again to record the disappearance of the Korotkoff sound as diastolic pressure.
- After measurement is completed, fully open the valve to release any remaining air in the cuff. The heartbeat rate is displayed on the LCD. (Fig. ⑤)

- Note:
 1. Sit upright in a chair, and take 5-6 deep breaths. Avoid leaning back while the measurement is being taken.
 2. Do not cross the legs while sitting and keep the feet flat on the floor during measurement.

Error Messages

- Er** Heartbeat rate detection error. The cuff is deflated too rapidly. Please keep the deflation rate among 2-3 mmHg per second.
EE A pressure of over 15 mmHg remaining for more than 3 minutes is detected. Please press the On/Off/Zeroing key to switch the cuff off, or the unit will automatically switch off after 60 seconds.
Data Error: Please remove and reload the batteries. If error keeps occurring, return the device to the local distributor or service center.
E3 The cuff is inflated to a pressure over 301 mmHg. Please open the valve immediately to rapidly exhaust the air in the cuff.

Maintenance

If the unit is used for over 10,000 times, the calibration reminding message (CAL) will appear every time the unit is switched on. It is recommended that calibration be completed every two years. Only the manufacturer or an authorized service technician is qualified for the calibration service. Please contact the manufacturer, Rossmax International Ltd. or your local distributors. (The PC Link socket for calibration purpose is inside the battery compartment.)

Safety Information

- Rossmax GA series are exclusively designed for blood pressure measurement at the upper arm or thigh on healthy skin.
- Please do not inflate to more than 300 mmHg.
- Do not take your blood pressure for more than 3 minutes.
- Wait for at least 5 minutes before repeated measurements.

Cautionary Notes

- The unit contains high-pressure assemblies. Therefore, avoid extreme temperatures, humidity, and direct sunlight. Avoid dropping or strongly shocking the main unit, and protect it from dust.
- Clear the blood pressure monitor body and the cuff carefully with a slightly damp, soft cloth. Do not press. Do not wash the cuff or use chemical cleaner on it. Never use thinner, alcohol or petrol (gasoline) as cleaner.
- Leaky batteries can damage the unit. Remove the batteries when the unit is not used for a long time.
- The unit should not be operated by children so to avoid hazardous situations.
- If the unit is stored near freezing, allow it to acclimate at room temperature before use.
- The unit is not field serviceable. You should not use any tool to open the device nor should you attempt to adjust anything inside the device. If you have any problems, please contact who you purchased GA series from or please contact Rossmax International Ltd.
- To stop operation at any time, open the air release valve, and the air in the cuff will be rapidly exhausted.
- Please note that this is a home healthcare product only and it is not intended to serve as a substitute for the advice of a physician or medical professional.
- Do not use this device for diagnosis or treatment of any health problem or disease. Measurement results are for reference only. Consult a healthcare professional for interpretation of pressure measurements. Contact your physician if you have or suspect any medical problem. Do not change your medications without the advice of your physician or healthcare professional.
- Electromagnetic interference. The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These may lead to temporary impairment of measurement accuracy.
- Dispose of device, batteries, components and accessories according to local regulations.
- This monitor may not meet its performance specification if stored or used outside temperature and humidity ranges specified in Specifications.
- Please note that when inflating, the functions of the limb in question may be impaired.
- Avoid any mechanical restriction, compression or bending of the cuff line.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous (AV) shunt.
- Do not apply the cuff on the side, where a mastectomy has been performed in your patient history.
- Do not place the cuff over wounds as this may cause further injury.
- Do not use the tubing and/or AC adaptor for any other purpose than those specified, as they can cause risk of strangulation.
- Do not service or maintain device and cuff while in use.
- During the blood pressure measurement, blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.
- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.

Specifications

Power Source	DC 3V Two AAA Batteries
Measurement Range	Pressure: 0-300 mmHg; Pulse: 30-240 beats/min
Accuracy	Pressure: 10% -40°C ±3 mmHg; Pulse: ±5% of reading
Pressure Sensor	Semi-conductor
Inflation	Manual Inflation
Deflation	Automatically operated deflation valves
Auto-shut-off	1 min. after last key operation
Back light auto-shut-off	15 sec. after last key operation
Operation environment	10°C -40°C (50°F -104°F); 85% RH max
Storage environment	10°C -40°C (50°F -104°F); 85% RH max
Dimensions	76(L) x 63(W) x 48(H) mm
Weight	90g (without batteries)
Limited users	Adult users

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
Specifications are subject to change without notice.

- GA101: without stethoscope
- GA102: with stethoscope
- GA112: with double head stethoscope

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