Nederland

Introductie

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accurate reading and avoid cross contamination approximately, no matter the probe cover has been installed or not.

- Do not drink hot or cold fluids, exercise, and smoke or perform other activities prior to a reading. These activities will raise or lower temperature readings when compared to your normal, average temperature.
- Please ensure the close contact between the probe and the armpit of the baby or children.
- Do not service or maintain device while in use.

Specifications

Measurement Range	32.0°C~42.9°C (90.0°F~109.9°F)
Accuracy	±0.1°C, 35.5°C~42.0°C (±0.2°F, 95.9°F~107.6°F) at standard room temperature 18°C~28°C (64.4°F~82.4°F); ±0.2°C at 35.5°C or over 42.0°C (±0.4°F over 95.9°F or over 107.6°F)
Display	Liquid crystal display, 3 1/2 digits
Battery	One 1.5V DC button size battery (SR41 or LR41)
Power consumption	0.15 milliwatts in measurement mode
Battery Life	Approx. 1500 times operation or 1 year with 1-2 measurement per day including standby mode.
Dimension	12.5cm x 2.1cm x 1.2cm (PxLxH)
Weight	Approx. 12.8 grams including battery
Alarm	Approx. 10 seconds sound signal when peak temperature reached
Environmental	Temperature: 5.0°C~40.0°C (41.0°F~104.0°F), Humidity: 15%~95%RH; 700~1060 hPa
Storage and transportation condition	Temperature: -20.0°C~55.0°C (-4.0°F~131.0°F), Humidity: 15%~95%RH; 700~1060 hPa
IP classification	IP22, Protection against harmful ingress of water and particulate matter
Safety classification	IEC 60601-1:2005+A1:2010, IEC 60601-2-56:2017+A1:2018 Medical electrical equipment - Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement.

Attention: Consult the accompanying documents. Please read this manual carefully before use. Please be sure to keep this manual.

Battery Replacement

Place a new 1.5V D.C. button size battery type SR41 or LR41 or equivalent in the chamber with positive side faced up and negative side faced down.

Display on the LCD

M: Last measured temperature.

N: When the appears in the bottom right corner of the LCD, the battery is exhausted and needs replacement.

Troubleshooting

Error message	Problem	Solution
	The system is not functioning properly	Unload the battery, wait for 1 min and repower it. If the message reappears, contact the retailer for service.
H -- C (bi--bi--bi)x10	Temperature taken is higher than 42.9°C (109.9°F)	Take a new temperature with close contact and sufficient test.
H -- C	Temperature taken is lower than 32.0°C (90.0°F)	Take a new temperature with close contact and sufficient test.
	Dead battery. Battery icon is flashing, can't be measurable.	Suggest to replace the battery.

Cleaning and Disinfection

1. Clean the unit by wiping it with a dry cloth and disinfect the probe with ethyl alcohol.

2. Don't let the unit contact any chemical thinner.

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• Do not store the unit under direct sunlight, at a high temperature, in high humidity or dust.

• Do not disassemble. See BATTERY REPLACEMENT to replace battery.

• Avoid strong electromagnetic interference such as microwave ovens and cell phones.

• Keep battery away from children. Harmful is swallowed.

• Battery should not be charged or placed into extreme heat as it may explode.

• Remove battery from the thermometer when not in operation for a long time.

• If the unit is stored in an extreme environment before measurement, please make sure it has been firstly adapted in the room temperature, so that the unit can be measured normally.

• Measurement results are for reference only. Contact your physician if you have or suspect any medical problems. Do not change your medications without the advice of your physician or healthcare professional.

• This device may not meet its performance specification if stored or used outside temperature and humidity ranges specified in specifications.

• If this device is used according to the operations instruction, periodic re-calibration is not required. If you still have questions, please send the device to dealers.

• Please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at your local retailer or at appropriate collection points provided in your country.

• High, prolonged fever requires medical attention especially for young children. Please contact your physician.

• For safe reason, during children's temperature measurement, please keep them from crying, walking, talking and any related dangerous activities.

• Usage of the probe may result in a 0.1°C(0.2°F) difference from actual temperature.

• To clean the probe before and after using the thermometer to ensure an

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