

- 1. Arm Cuff
- 2. LCD Display
- 3. Air Tube and Connector
- 4. ON/OFF/START Key
- 5. Battery Cover
- 6. Data Link Socket
- 7. AC Adaptor Jack

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EN English

Introduction
Blood pressure measurements determined with X1 are equivalent to those obtained by a trained operator using a cuff/stethoscope auscultation method, when the limits prescribed by the American National Standard, Electronic or Automated Sphygmomanometry. This unit is to be used by adult consumers in a home environment. The patient is an intended operator. Do not use this device on infants or neonates. X1 is protected against manufacturing defects by an established International Warranty Program. For warranty information, you can contact the manufacturer, Rossmax International Ltd. Attention: Consult the accompanying documents. Please read this manual carefully before use. For specific information on your own blood pressure, contact your physician. Please be sure to keep this manual.

Real Time Measuring Technology
This unit uses the oscillometric method to measure your blood pressure. Before the cuff starts inflating, the device will establish a baseline cuff pressure equivalent to the air pressure. This unit will determine the appropriate inflation level based on pressure oscillations, followed by cuff deflation. During the deflation, the device will detect the amplitude and slope of the pressure oscillations and thereby determine for you the systolic blood pressure, diastolic blood pressure, and pulse.

Preliminary Remarks
This Blood Pressure Monitor complies with the European regulations and bears the CE mark, "CE 1639". The quality of the device has been verified and conforms to the provisions of the Council Directive 93/42/EEC (Medical Device Directive), Annex I essential regulated and applied harmonized standards.

EN 1060-1: 1995/A2: 2009 Non-invasive sphygmomanometers - Part 1 - General requirements.
EN 1060-3: 1997/A2: 2009 Non-invasive sphygmomanometers - Part 3 - Supplementary requirements for the design and construction of the device.
EN 1060-4: 2004 Non-invasive sphygmomanometers - Part 4 - Test Procedures to determine the in-vitro in-vivo accuracy of automated non-invasive sphygmomanometers. ISO 81060-2:2013 In-vivo accuracy of automated non-invasive sphygmomanometers. Part 2: Clinical investigation of automated measurement type.

The blood pressure monitor is designed for long service time. Ensure continued accuracy. It's recommended that digital blood pressure monitors require recalibration. This monitor (under normal usage with approx. 3 measurements a day) does not require recalibration for 2 years. Once the unit should be recalibrated the device will display "CR". The unit should also be recalibrated if the monitor sustains damage due to bursted cuff or inflating tube and/or extreme hot or cold temperature/humidity changes. When "CR" appears, simply return to your nearest dealer for recalibration service.

Blood Pressure Standard
Refer to the definitions of the World Health Organization, the blood pressure ranges can be classified as follows: (Ref. 1999 WHO International Hypertension Guidelines for the management of hypertension). This blood pressure classification are based on historical data, and may not be directly applicable to any particular patient. It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk. For reliable monitoring of blood pressure, keeping log records is recommended. Please download the blood pressure log at our website www.rossmax.com.

Display Explanations
EE / Measurement Error: Make sure the L-plug is securely connected to the air socket and the cuff correctly and keep arm steady during measurement. If the error keeps occurring, return the device to your local distributor or service center.
E1 / Air Circ Abnormality: Make sure the L-plug is securely connected to the air socket on the side of the unit and measure again quietly. If the errors still occur, return the device to your local distributor or service center for help.
E2 / Pressure Exceeding 300 mmHg: Switch the unit to the error measure again quietly, if the error keeps occurring, return the device to your local distributor or service center.
E3 / Data Error: Remove the batteries, wait for 60 seconds, and reload. If the error occurs again, return the device to your local distributor or service center.
Er / Exceeding Measurement Range: Measure again quietly, if the error keeps occurring, return the device to your local distributor or service center.

Cuff Wrap Detection
If the cuff was wrapped too loosely, it may cause unreliable measurement results. The "Cuff" icon will help you to determine if the cuff is wrapped snugly enough. The specified icon appears once a "loosen cuff" has been detected during measurement. Otherwise the specified icon appears if the cuff is wrapped correctly during measurement.

Movement Detection
The "Movement Detection" helps reminding the user to remain still and is indicating movement during measurement. The specified icon appears once a "body movement" has been detected during and after each measurement.

Hypertension Risk Indication (HRI)
The World Health Organization, classifying blood pressure ranges into 6 grades. This unit is equipped with innovative blood pressure risk indication, which visually indicates the risk level (optimal / high-normal / grade1 hypertension / grade2 hypertension / grade3 hypertension) of the result after each measurement.

Irregular Heartbeat (IHB) Detection
This unit is equipped with an Irregular Heartbeat (IHB) Detection which allows those who have an irregular heartbeat to obtain accurate measurements alerting the user of the irregular heartbeat during the measurement.
Note: It is strongly recommended that you consult your physician if the IHB icon appears often.

Using the AC Adaptor (Optional)
1. Connect the AC adaptor with the AC adaptor jack on the right side of the unit.
2. Connect the AC adaptor with the power cord. (AC adaptors with required voltage and current indicated near the AC adaptor jack.)
3. Please unplug the batteries when operating with the AC mode for a longer period of time. Leaving the batteries in the compartment for a long time may cause leakage, which may harm to the batteries of the unit.
4. 2 batteries are needed when operating with the AC mode.
5. AC adaptors are optional. Please contact the distributor for the compatible AC adaptors.
6. Use only the authorized AC Adaptor for this blood pressure monitor. Information for the authorized AC adaptor, please refer to APPENDIX 1.
7. Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.
8. Do not use the tubing and/or AC adaptor for any other purpose than those specified, as they can cause risk of strangulation.
9. Do not service or maintain device and cuff while in use.

Measurement Method
Measurement Range: 0-260 mmHg; Pulse: 40-199 beats/minute
Pressure Sensor: Semi conductor
Accuracy: Precision: ± 3 mmHg; Pulse: ± 5% of reading
Pump Drive: Automatic Air Release Valve
Memory capacity: Last Number Memory Recall
Auto-shut-off: 1 minute after last key operation
Permissible Operating Temp: 10°C-40°C (50°F-104°F); 15%-85% RH
Relative and Humidity: 10%-90% RH
Permissible Transport and Store: -10°C-60°C (14°F-140°F); 10%-90% RH; age Temperature and Humidity 700-1060 hPa

Applying the Cuff
1. Unwrap the arm cuff, leaving the end of the cuff through the D-ring of the cuff.
2. Put your left arm through the cuff. The color strip indication should be positioned closer to you with the tube pointing in the direction of your arm (Fig. ①). Turn your left palm upward and place the edge of the arm cuff at approximately 1.5 to 2.5 cm above the inner side of the elbow joint (Fig. ②). Tighten the cuff by pulling the end of the cuff.
3. Center the tube over the middle of the arm. Press the hook and loop material together securely. Allow room for 2 fingers to fit between the cuff and your arm. Position the artery mark (③) over the main artery (on the inside of your arm) (Fig. ③, ④). Note: Locate the main artery by pressing with 2 fingers approximately 2 cm above the bend of your elbow on the upper left arm. Identify where the pulse can be felt the strongest. This is your main artery.
4. Plug the cuff connecting tube into the cuff (Fig. ⑤).
5. If your arm or table (palm) surface and the cuff is at the same height as your heart. Make sure the tube is not kinked (Fig. ⑥).
6. This cuff is suitable for you if the arrow falls within the solid color line as shown on the right (Fig. ⑦). If the arrow falls outside the solid color line, you will need a cuff with other circumferences. Contact your local dealer for additional size cuffs.

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The blood pressure monitor is designed for long service time. Ensure continued accuracy. It's recommended that digital blood pressure monitors require recalibration. This monitor (under normal usage with approx. 3 measurements a day) does not require recalibration for 2 years. Once the unit should be recalibrated the device will display "CR". The unit should also be recalibrated if the monitor sustains damage due to bursted cuff or inflating tube and/or extreme hot or cold temperature/humidity changes. When "CR" appears, simply return to your nearest dealer for recalibration service.
Blood Pressure Standard
Refer to the definitions of the World Health Organization, the blood pressure ranges can be classified as follows: (Ref. 1999 WHO International Hypertension Guidelines for the management of hypertension). This blood pressure classification are based on historical data, and may not be directly applicable to any particular patient. It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk. For reliable monitoring of blood pressure, keeping log records is recommended. Please download the blood pressure log at our website www.rossmax.com.

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Cuff Wrap Detection
If the cuff was wrapped too loosely, it may cause unreliable measurement results. The "Cuff" icon will help you to determine if the cuff is wrapped snugly enough. The specified icon appears once a "loosen cuff" has been detected during measurement. Otherwise the specified icon appears if the cuff is wrapped correctly during measurement.
Movement Detection
The "Movement Detection" helps reminding the user to remain still and is indicating movement during measurement. The specified icon appears once a "body movement" has been detected during and after each measurement.
Hypertension Risk Indication (HRI)
The World Health Organization, classifying blood pressure ranges into 6 grades. This unit is equipped with innovative blood pressure risk indication, which visually indicates the risk level (optimal / high-normal / grade1 hypertension / grade2 hypertension / grade3 hypertension) of the result after each measurement.
Irregular Heartbeat (IHB) Detection
This unit is equipped with an Irregular Heartbeat (IHB) Detection which allows those who have an irregular heartbeat to obtain accurate measurements alerting the user of the irregular heartbeat during the measurement.
Note: It is strongly recommended that you consult your physician if the IHB icon appears often.

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