



"DK" Deluxe Automatic Blood Pressure Monitor

Model: AC701 **k**



- Digital Korotkoff Technology
- Hypertension Risk Indication
- Irregular Heartbeat Detection (IHB)
- 120 memories for 2-users with date and time indication
- Average of last 3 readings
- Jumbo size "XL" digit
- One-touch automatic operation
- Carrying bag
- Patented Digital Korotkoff Cone cuff
- Lifetime calibration
- PC Link (Optional)
- Adaptor socket



Digital Korotkoff sound (DK) Technology



Hypertension Risk Indication



Irregular Heartbeat Detection (IHB)



120 memories for 2-users with date and time indication



Average of last 3 readings




Jumbo size "XL" digit



AAMI

This blood pressure monitor meets the AAMI accuracy standard requirements (AAMI/ANSI/ISO 81060-2)



Cone cuff fits arm circumference: 24~36cm (9.4"~14.2")
Other size cuffs available upon request

«DK»® Digital Korotkoff Technology

With a proprietary microphone built into the arm cuff, the innovative DK Technology actually listens like your doctor when taking a reading. The Exclusive **DK Technology** delivers **ABSOLUTE Accuracy** by using the advanced dual sensor cuff system that senses both pressure and sound.

Model	Qty per carton	Carton volume
AC701k	8 pcs	0.058 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
9882 pcs	19764 pcs	22932 pcs

○ Storage and Transportation Condition
-10°C~60°C(14°F~140°F)
10%~90% RH, 700~1060 hPa