



## Air Mattress Bubble

Model: AM30



Adjustable hangers to suit any kind of hospital bed frame

2 to 1 alternating therapy


- Bubble mattress
- 2 to 1 alternating therapy
- Adjustable hangers for most bed types
- Durable medical grade PVC
- Extension flags for stabilization
- High quality pump
- Weight indicator pressure adjustment setting
- 6 minutes cycling time
- Near silent low noise
- For prevention and treatment of patients at low risk of pressure ulcers



Low noise design




Medical grade PVC



2 to 1 alternating therapy



Extension flags for stabilization



Easy to use adjustment



Adjustable hangers to suit any kind of hospital bed frame

### Accessories



High Quality Rubber Tube



Fuse





Repairing Kit



Mattress

**REF**

Fuse: AM\_AC\_006\_00, Rubber tube: AM\_AC\_004\_00  
Mattress: AM\_MA\_003\_00, Repairing kit: AM\_AC\_007\_00



Model	Qty per carton	Carton volume
AM 30	4 pcs	0.0408 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
2744 pcs	5492 pcs	6372 pcs

Storage and Transportation Condition  
-18°C to 43°C (0°F to 110°F)  
10~90% RH

www.rossmax.com

Epaper\_AM30\_ver2122