

## Air Mattress Cell

Model: AM40



Adjustable hangers to suit any kind of hospital bed frame

2 to 1 alternating therapy


- Cell mattress
- 2 to 1 alternating therapy
- Adjustable hangers for most bed types
- Comfortable Nylon mattress
- High quality pump
- Weight indicator pressure adjustment setting
- 12 minutes cycling time
- Near silent low noise
- For treatment of patients at medium to high risk of pressure ulcers




Low noise design



Comfortable Nylon mattress



2 to 1 alternating therapy



Easy to use adjustment



Adjustable hangers to suit any kind of hospital bed frame

### Accessories



High Quality Rubber Tube



Spare Cell




Fuse




Cover

REF

Fuse: AM\_AC\_005\_00, Rubber tube: AM\_AC\_003\_00  
Cell: AM\_MA\_002\_00, Cover: AM\_AC\_002\_00



● A ● B



Model	Qty per carton	Carton volume
AM 40	2 pcs	0.0579 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
1000 pcs	2000 pcs	2246 pcs

Storage and Transportation Condition  
-18°C to 43°C (0°F to 110°F)  
10~90% RH

www.rossmax.com

Epaper\_AM40\_ver2212