



Aero Spacer (Valved Holding Chamber)

Model: AS175



AS175 with Cap



AS175 with L-size Mask



AS175 with M-size Mask



AS175 with S-size Mask

- High efficiency one way valve for easy inhalation and exhalation
- Universal Backpiece, can be used with all commonly prescribed Metered Dose Inhalers
- Anti-static spacer improves the delivery of MDI by inhibiting the adherence of drug to spacer wall
- Soft and flexible silicone mask provides comfortable touch and reliable seal
- Portable design & easy holding
- Easy disassembled and reassembled for cleaning
- Mouthpiece cap can protect the mouthpiece from dirt and debris when the device is not in use
- This product contains no latex and BPA free.
- Small, medium & large facemask to cater to the different needs
- Small facemask: 0-18 months
- Medium facemask: 12 months – 5 yrs.
- Large facemask or mouthpiece: 5 yrs. +
- Chamber volume: 175ml
- For single patient use only

AS175	Qty/ctn	CBM/ctn	20FCL	40FCL	40HQ
With Mouthpiece & Cap (5 years +)	50	0.064	21,000	42,000	49,000
With S-size Mask (0-1.5 years)	50	0.064	21,000	42,000	49,000
With M-size Mask (1 – 5 years)	50	0.064	21,000	42,000	49,000
With L-size Mask (5 years +)	50	0.075	16,800	32,400	41,250

● Storage and Transportation Condition
5°C to -40°C (41°F to 104°F)
15-95% RH

www.rossmax.com

OBM_Epaper_AS175_ver1605