

AC Adaptor for Blood Pressure Monitor / Nebulizer

• US version



• EU version



Energy Saving



Regulated Output

For Blood Pressure Monitor


Model No.	HK-X205-A06 / HK-X105-A06
Input	100-240V, 50-60Hz
Output	DC 6V 0.8A
Plug size	Outer \ominus : \varnothing 4.0 mm Inner \oplus : \varnothing 1.7mm

Model No.	HK-X111-A075 / HK-X211-A075
Input	100-240V, 50-60Hz
Output	DC 7.5V 1.5A
Plug size	Outer \ominus : \varnothing 3.8 mm Inner \oplus : \varnothing 1.35 mm

For Nebulizer

Model No.	HK-X112-A12 / HK-X212-A12
Input	100-240V, 50-60Hz
Output	DC 12V 1.0A
Plug size	Outer \ominus : \varnothing 3.8 mm Inner \oplus : \varnothing 1.35 mm

Model No.	HK-X205-A06 / HK-X105-A06
Input	100-240V, 50-60Hz
Output	DC 6V 0.8A
Plug size	Outer \ominus : \varnothing 3.5 mm Inner \oplus : \varnothing 1.3 mm



Model	Qty per carton	Carton volume
Adaptor / US version	60 pcs	0.021 cbm/ctn
Adaptor / EU version	60 pcs	0.021 cbm/ctn

- Storage and Transportation Condition :
-20°C~60°C, \leq 90% RH, 700-1060 hPa
- Operating condition : 10°C~40°C,
 \leq 90% RH, 700-1060 hPa