





"XL" Deluxe Automatic Wrist Blood Pressure Monitor

Model: BQ705





- Jumbo size "XL" digit
- 120 memories for 2-users with date and time indication
- Guest Mode
- Movement Detection
- Smart sense technology
- Hypertension Risk Indication
- Irregular Heartbeat Detection (IHB)
- Average of last 3 readings
- One-touch automatic operation
- Carrying case
- PC Link (Optional)
- Lifetime calibration









Jumbo size "XL" digit

120 memories for 2-users

Guest Mode

Movement Detection









Smart Sense Technology

Hypertension Risk Indication

Irregular Heartbeat Detectjion (IHB)

Average of last 3 readings





This blood pressure monitor meets the AAMI accuracy standard requirements (AAMI/ANSI/ISO 81060-2)

Smart Sense Technology

Smart Sense Technology utilizes an advanced air control system. The electronic deflation valve precisely controls the deflation rate and allows the monitor to detect and analyze pulse signals, adjusting cuff pressure during deflation. This significantly ensures more precise analysis, while reaching greater accuracy.



Cuff fits wrist circumference: 13.5~22cm (5.3"~8.7")

Model	Qty per carton	Carton volume
BQ705	36 pcs	0.046 cbm/ctn

Storage and Transportation Condition -10°C~60°C(14°F~140°F) 10%~90% RH, 700~1060 hPa