



Digital Sphygmomanometer

Model : **GA series**
Manduas II



Rossmax patented Digital Sphygmomanometer. Always maintain calibration even if dropped. Spill proof design with backlight. 2 AAA batteries provide up to 3,000 ements. Built-in alert feature allows easy to program calibration after 10,000 measurements.

- ESH validated for professional clinical use
- Extra large backlight LCD
- Shock and splash resistance
- Pulse rate display
- Virtually indestructible design
- 360° swivel cuff clip
- Patent No.: CN: ZL01259708.2; TW: 488518; GE: 20206337.2; JP: 3089650; US: 717746



Marker



Extra large backlight LCD



Shock and splash resistance



360° swivel cuff clip



Chrome-plate brass air release valve

GA102 - D-ring cuff with single head stethoscope



Cuff fits arm circumference:
24~36cm (9.4" ~14.2")

Model	Qty per carton	Carton volume
GA series	18 pcs	0.067 cbm/ctn

Storage and Transportation Condition
-34°C ~60°C (30°F ~149°F)
10%~90% RH, 700~1060 hPa