



rossmax

just a heartbeat away

prevention • monitoring • therapy • professional

Aneroid Sphygmomanometer

Model : **GB series**
AGC



GB100



- White enamel 300 mmHg no-pin stop manometer
- Chrome-plate brass air release valve
- Latex free Cuff
- Single head Stethoscope



White enamel 300mmHg non-stop pin reliable manometer



Chrome-plate brass air release valve



Light-weight stethoscope



Latex free cuff

GB102 - D-ring cuff with single head stethoscope



Rossmax Cuff with D-ring :

Adult: 24 - 32cm (9.4" - 12.6")

Other size cuffs available upon request

Model	Qty per carton	Carton volume
GB series	20 pcs	0.042 cbm/ctn

- Storage and Transportation Condition
-10°C~60°C(14°F~140°F)
10%~90% RH, 700~1060 hPa

www.rossmax.com
OBM_Epaper_GBseries_ver1546