





Infrared Non-Contact Temple Thermometer

Model: HA500



- Measurement distance within 5cm (2 inch)
- 1 second measurement
- Temple Temperature
- Object temperature
- Fever alarm
- 9 memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction



Patent Nos : ZL 2011 3 0359484.6

Clinically Accurate

ASTM E1965-98:2009

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.



Model	Qty per carton	Carton volume
HA500	50 pcs	0.072 cbm/ctn

Storage and Transportation Condition
-25°C~55°C(-13°F~131°F)
≤95% RH, 700~1060 hPa

www.rossmax.com
OBM_Epaper_HA500_ver1546