



Infrared Non-Contact Temple Thermometer

Model: HA500



Temple Thermometer Non Contact



- Measurement distance within 5cm
- 1 second measurement
- Temple Temperature
- Object temperature
- Fever alarm
- 9 memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction



Patent Nos : ZL 2011 3 0359484.6

Clinically Accurate

ASTM E1965-98:2009

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

| Model | Qty per carton | Carton volume |
|-------|----------------|---------------|
| HA500 | 50 pcs | 0.072 cbm/ctn |

● Storage and Transportation Condition
-25°C~55°C(-13°F~131°F)
≤95% RH, 700~1060 hPa