




Infrared Non-Contact Temple Thermometer

Model: HC700



- Telephoto temperature measurement distance within 10cm (3.94 inch)
- 1 second measurement
- Temple temperature
- Object temperature
- Auditory fever alarm
- 9 memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction



Patent Nos :
ZL 2014 3 0054262.7

Clinically Accurate

ASTM E1965-98:2009

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

Model	Qty per carton	Carton volume
HC700	30 pcs	0.060 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
14,010 pcs	28,020 pcs	32,520 pcs

- Storage and Transportation Condition
-25°C ~ 55°C (-13°F ~ 131°F)
≤ 85% RH, 700~1060 hPa