



## Deluxe Automatic Blood Pressure Monitor

Model: MJ701

### Movement Detection



The "Movement Detection" helps reminding the user to remain still and is indicating any body movement during measurement. The specified icon appears once a "body movement" has been detected during and after each measurement.



- Movement Detection
- Real Fuzzy Technology
- Hypertension Risk Indication
- Irregular Heartbeat Detection (IHB)
- Average of last 3 readings
- 120 memories for 2-users with date and time indication
- Guest mode
- Latex-free patented universal cone cuff
- PC Link (Optional)
- Adaptor socket
- Lifetime Calibration
- PC Link (Optional)
- Carrying bag



Movement Detection



Real Fuzzy Technology



Hypertension Risk Indication



Irregular Heartbeat Detection (IHB)



Average of last 3 readings



120 memories for 2-Users with Date and Time Indication



Guest Mode



Universal Cone Cuff



Dual Power Ready



Lifetime Calibration



### Clinically Validated

This blood pressure monitor has passed the criteria of the European Society of Hypertension by achieving all requirements and achieved the highest possible grade (A/A) according to the criteria of the British Hypertension Society. This monitor further complies with the standard required by the Association for the Advancement of Medical Instrumentation (AAMI).



Universal Cone cuff fits arm circumference: 24-40cm (9.4" ~ 15.7")  
Other size cuffs available upon request

Patent Nos:  
US: 7597667B2  
TW: 191076; TW: 1 236891; TW: I322262  
CN: 469964, Other patents pending

### Real uzzy Technology

Real Fuzzy Technology determines ideal cuff pressure based on one's systolic blood pressure and arm size. Users no longer need to pre-set the inflation level before measurement. The technology eliminates arm discomfort through unnecessary high inflation settings and also prevents inaccurate readings caused by erroneous cuff inflation levels.

| Model       | Qty per carton | Carton volume |
|-------------|----------------|---------------|
| MJ701f      | 12 pcs         | 0.047cbm/ctn  |
| QTY per 20' | QTY per 40'    | QTY per 40'HQ |
| 7116 pcs    | 14232 pcs      | 16524 pcs     |

- Storage and Transportation Condition  
-10°C ~ 60°C (14°F ~ 140°F)  
10% ~ 90% RH, 700 ~ 1060 hPa