


## Nebulizer Kit\_2in1 Cup

Model: N4





• 2in1 Cup with Mouthpiece

- 2in1 cup for easy clean, assembly and disassembly
- Your nebulizer kit is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy.
- Single patient use
- Nebulizer kit should be replaced after using for 12 months.
- Storage temperature range: -4°F to 140°F (-20°C to 60°C)
- Consistent particle size (MMAD) less than 3µm
- Average nebulization rate



### Hygienic 2in1 Cup Design



Model	Qty per carton	Carton volume
Nebulizer kit with Mouthpiece	160 pcs	0.072 cbm/ctn

● Storage and Transportation Condition  
-20°C~60°C(-4°F~140°F)  
≤ 90% RH; 700~1060 hPa

www.rossmax.com

OBM\_Epaper\_N4\_ver1516