



Infrared Ear Thermometer

Model : RA600 Cutie



Storage Case



Patent Nos : ZL 2011 3 0359483.1
ZL 2010 2 0668110.2
M403995
Other patents pending



Clinically Accurate

This Ear thermometer has been clinically tested in several big teaching hospitals based on the protocol of EN 12470-5:2003 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

EN 12470-5: 2003

- Highest accuracy by rotating in 3" Scanning Mode or 1" measurement
- Multiple- Use PLUG & SECURE Probe Cover
- Fever alarm
- 9 memories
- C/F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction

RESULT IN 1" OR 3" Scanning Mode



Model	Qty per carton	Carton volume
RA600 Cutie	24 pcs	0.044 cbm/ctn

Storage and Transportation Condition
-25°C~55°C(-13°F~131°F)
≤95% RH, 700~1060 hPa

www.rossmax.com

OBM_Epaper_RA600Q_ver1907