

Fingertip Pulse Oximeter

Model: SB100




- %SpO₂ & Heart Rate
Big LED screen display for easy reading
- Shield Design against ambient light
- Biocompatibility protection
- One-touch ON/OFF key
- Pulse Strength Indicator
- String-attached
Economical & portable health management tool



%SpO₂ & Heart Rate



Shield Design against Ambient Light



Biocompatibility protection



Pulse Strength Indication



String-Attached



Large Font Size

Specification

- Measurement range: SpO₂: 35%~99%; HR: 30~250 bpm
- Precision: SpO₂: 70%~99%; ±2%; 35%~69% (unspecified); Pulse rate: 30~250bpm; ±3
- Power Supply: 2"AAA" Alkaline

Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring

Model	Qty per carton	Carton volume
SB100	40 pcs	0.028 cbm/ctn

- Storage temperature: -25°C - 70°C (-13°F - 158°F)
- Relative humidity: 15% - 90% (no condensing)
- Atmospheric Pressure: 700 to 1013 hPa