

Thermo Cal



Function

- Check the accuracy of thermometer

Specification

Power range	AC100~120V 50/60 Hz AC210~230V 50/60 Hz	
Power of water heater	255W	
Testing temperature	Ambient temperature+ 5°C to maximum 45°C	
Ambient temperature	16 ~ 30°C	
Stability of water bath	± 0.02°C	
Diameter of blackbody cavity	Pen type	15 mm
	Ear type	12.5 mm
	Forehead type	55 mm
Capacity	5.39 Liter	
Weight	9.69kgs	
Dimension	180(W)x350(H)x320(D) mm	



Pen Type



Ear Type



Forehead Type

Rossmax Thermo Cal

Accuracy is always the biggest concern of a thermometer user. Rossmax' Thermo Cal is a portable and convenient device that allows the user to easily check the accuracy of a thermometer in just a few seconds. The innovative design of the adjustable blackbody cavity applies to all three types of Rossmax thermometer products - pen type, ear and forehead thermometer.

Model	Qty per wooden carton	Carton volume
Thermo Cal	1 pcs	0.115 cbm/ctn