



Comfort Suction Unit

Model: V3

Triple Protection consist of

1 Overflow valve, 2 Antibacterial filter 3 Fool-proof design



Patent Nos:
ZL 2016 3 0298084.1
Other patents pending

It's Comfort

- High performance pump delivers -550mmHg vacuum and 16 liter/min flow rate

It's Protected "Triple Protection"

- Overflow valve plus antibacterial filter and fool-proof design that prevents back flow of fluid and airborne contaminants from entering the pump

It's Convenient

- Easy to store using tubing holder and accessory compartment

It's Ergonomic

- Comprehensive dashboard with ergonomic view angle ensures a smooth operation

Technical Specification

- Suction class: High Vacuum / Low Flow
- Max vacuum @ pump inlet: -73 kPa/ -550 mmHg
- Max air flow @ pump inlet: 16 L/min
- Adaptor Input: AC 100-240V, 50/60Hz
- Adaptor Output: DC 12V 2.0A
- Collection canister capacity: 1400ml (Polycarbonate, PC)
- Power Consumption: ≤ 24W
- * Subject to technical modification without prior notice

Tubing holder



Accessory compartment



Ergonomic design



Model	Qty per carton	Carton volume
V3	4 pcs	0.131 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
480 pcs	960 pcs	960 pcs
QTY per pallet	48 pcs	