


## Baby Scale

Model : WE300 Cutie



- Easy to read
- Tare function
- Discretional hold function
- Auto shut-off
- Low battery indicator
- Overload indicator
- Unit switchable: kg / lb / oz
- Maximum capacity: 20kg
- Minimum capacity: 100g
- Graduation: 5g



Model	Q'ty per carton	Carton volume
WE300	4 pcs	0.0763 cbm/ctn
Q'ty per 20'	Q'ty per 40'	Q'ty per 40' HQ
1468 pcs	2936 pcs	3408 pcs