




## Premium Accessories Pack for HS200 - Blood Glucose Monitoring System



- 10 Test strips
- 10 Lancets
- Lancing device
- Control solution
- Zipper bag
- Instruction Manual
- Quick guide

### Advanced GDH-FAD Technology

The GDH-FAD test strips can eliminate interference from oxygen variations and non-glucose sugar such as maltose and galactose, providing you with the most reliable and accurate results.



**GDH-FAD**

Enzyme	GDH-FAD	GDH-PQQ	GOD
Affected by Oxygen	No	No	Yes
Affected by Galactose	No	Yes	No
Affected by Maltose	No	Yes	No

● Storage and Transportation Condition  
4°C ~ 30°C (39°F ~ 86°F)  
10% ~ 90% RH, 700 ~ 1060 hPa