



Test Strips for HS200 - Blood Glucose Monitoring System




- Advanced GDH- FAD Technology (applicable for wide HCT range: 0-70%)
- Measures with small Blood Sample
- 25 pieces or 50 pieces in a Box
- Not affected by Oxygen
- Not affected by Galactose
- Not affected by Maltose



Advanced GDH-FAD Technology

The GDH-FAD test strips can eliminate interference from oxygen variations and non-glucose sugar such as maltose and galactose, providing you with the most reliable and accurate results.



GDH-FAD

Enzyme	GDH-FAD	GDH-PQQ	GOD
Affected by Oxygen	No	No	Yes
Affected by Galactose	No	Yes	No
Affected by Maltose	No	Yes	No

● Storage and Transportation Condition
4°C ~ 30°C (39°F ~ 86°F)
10% ~ 90% RH, 700 ~ 1060 hPa