

Nebulizer Kit

Model: N2




Nebulizer Kit
Consisting of:
• N2 Bottle
• Mouthpiece

Specifications

- Consistent particle size (MMAD) less than $3.8\mu\text{m}$ (for NL100)
- Average nebulization rate $\geq 0.20\text{ml/min}$ (for NL100)

Use

- The nebulizer kit is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy.
- Single patient use
- The nebulizer kit should be replaced after using for 6 months.



Model	Qty per carton	Carton volume
Nebulizer Kit	160 pcs	0.072 cbm/ctn

● Storage and Transportation Condition
-20°C~60°C(-4°F~140°F)
≤ 90% RH; 700~1060 hPa

www.rossmax.com

OBM_Epaper_N2_ver2049