

## Nebulizer Kit

Model: N2




Nebulizer Kit  
Consisting of:  
• N2 Bottle  
• Mouthpiece

### Specifications

- Consistent particle size (MMAD) less than  $3.8\mu\text{m}$  (for NL100)
- Average nebulization rate  $\geq 0.20\text{ml/min}$  (for NL100)

### Use

- The nebulizer kit is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy.
- Single patient use
- The nebulizer kit should be replaced after using for 6 months.



Model	Qty per carton	Carton volume
Nebulizer Kit	160 pcs	0.072 cbm/ctn

● Storage and Transportation Condition  
-20°C~60°C(-4°F~140°F)  
≤ 90% RH; 700~1060 hPa

[www.rossmax.com](http://www.rossmax.com)

OBM\_Epaper\_N2\_ver2326