

Innovative Triggered Nebulizer Bottle

Model: N3



N3 Patent No. :
EP238035
D140311
M388940
M411285
ZL201020179254.1
ZL201120021244X

Nebulizer Kit
Consisting of:
• Trigggle Bottle (model N3)
• Mouthpiece

Innovative Triggered Nebulizer Bottle

Rossmax Triggered Nebulizer Bottle (N3) is an innovative two-in-one nebulizer that produces aerosol at user's discretion without taking the continuing-mode aerosol mechanism away from users. N3 helps to have the prescribed medication delivered deep into patient's lungs and stays as patient paced it; most important of all, without medication waste and to achieve a time-saving treatment.

Trigggle Bottle

Specifications

- Lock on - continuing mode
Lock off - discretionary operation
- Leak free
- Breath enhance technology (one-way valve)
- Storage temperature range: -4°F to 140°F (-20°C to 60°C)
- Consistent particle size (MMAD) less than 3 μm
- Ergonomic bottle design

Use


- The nebulizer kit is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy.
- Single patient use



Lock on -
continuing mode



Lock off -
discretionary operation



Model	Qty per carton	Carton volume
Trigggle Bottle	160 pcs	0.072 cbm/ctn

● Storage and Transportation Condition
-20°C~60°C(-4°F~140°F)
≤90% RH; 700~1060 hPa