

## Nebulizer Kit


Model: N4



N4 Patent Nos :  
ZL 2015 3 0231205.6  
Other patents pending

Nebulizer Kit  
Consisting of:  
• 2in1 cup (model N4)  
• Mouthpiece

### Hygienic "2in1 Cup" Design




### Hygienic "2in1 Cup"

#### Specifications

- 2in1 cup for easy cleaning, assembly and disassembly
- Storage temperature range: -4°F to 140°F (-20°C to 60°C)
- Consistent particle size (MMAD) less than 2.6µm
- Average nebulization rate ≥ 0.35ml/min

#### Use

- The nebulizer kit is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy.
- Single patient use
- The nebulizer kit should be replaced after using for 6 months.



Model	Qty per carton	Carton volume
Nebulizer Kit	160 pcs	0.072 cbm/ctn

● Storage and Transportation Condition  
-20°C~60°C(-4°F~140°F)  
≤ 90% RH; 700~1060 hPa

www.rossmax.com

OBM\_Epaper\_N4\_ver2049