

Nebulizer Kit

Model: N4



N4 Patent Nos :
ZL 2015 3 0231205.6
Other patents pending

Nebulizer Kit
Consisting of:
• 2in1 cup (model N4)
• Mouthpiece


Hygienic "2in1 Cup"

Specifications



- 2in1 cup for easy cleaning, assembly and disassembly
- Storage temperature range: -4°F to 140°F (-20°C to 60°C)
- Consistent particle size (MMAD) less than 2.6µm
- Average nebulization rate ≥ 0.35ml/min

Use

- The nebulizer kit is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy.
- Single patient use
- The nebulizer kit should be replaced after using for 6 months.



Hygienic "2in1 Cup" Design



Model	Qty per carton	Carton volume
Nebulizer Kit	160 pcs	0.072 cbm/ctn

● Storage and Transportation Condition
-20°C~-60°C(-4°F~140°F)
≤ 90% RH; 700~1060 hPa

www.rossmax.com

OBM_Epaper_N4_ver1714