

## Suction Tester



### Suction Tester


Rossmax Suction Tester is a portable and easy-to-use device that allows the users to quickly check the performance of the pump of a suction unit. Power source is not required for Suction Tester during operation and the test result is easy to read and interpret.

This device is used to check the working conditions of the pump.


### Check

1. Maximum (free) air flow
2. Maximum vacuum

### Suction tester with Suction unit



### Suction tester with NK1000 (Nasal Aspirator)



### Specification

<b>Power source</b>		No need
<b>Operating temperature</b>		0 ~ 54°C
<b>Storage temperature</b>		-10 ~ 60°C
<b>Air flow</b>	scale	0~50 L/min
	accuracy	±4% of full scale
<b>Vacuum</b>	scale	0~76 cmHg
	accuracy	±2.5% of full scale
<b>Weight</b>		730g
<b>Dimension</b>		120(L) x 105(W) x 145(H) mm

Model	Qty per carton	Carton volume
Suction Tester	1 pcs	0.007 cbm/ctn