



rossmax

YOUR TOTAL HEALTHSTYLE PROVIDER



www.rossmax.com

YOUR TOTAL HEALTHSTYLE PROVIDER
TECHNOLOGIES
PORTFOLIO
MARKETING
AFTER-SALES
VALIDATION
CONTACT US

www.rossmax.com

In today's environmentally conscious world, technology is rapidly changing, not only in the health industry but also how we are delivering the right information to the right people at the right time. And in today's fast-paced world, paper just can't keep up! That's why Rossmax decided to delivering vital business information in a fast,



prevention



HEALTHSTYLE APP
SCALE

monitoring



THERMOMETER
BLOOD PRESSURE MONITOR
PULSE OXIMETER
BLOOD GLUCOSE METER

therapy



NEBULIZER
ASPIRATOR
PEAK FLOW METER & AERO SPACER
HEATING PAD

professional



HYPERTENSION MANAGEMENT
DIABETES MANAGEMENT
FEVER MANAGEMENT
RESPIRATORY CARE
WOUND HEALING

effective and highly efficient manner - namely through www.rossmax.com. Through monthly newsletters and social media platforms we communicate with friends and followers daily. Join the Rossmax digital world and at the same time help to protect the environment!

Rossmax' portfolio of advanced technologies, products and services provide health improvement for everyone.

Quality solutions for babies, kids, adults & elderly people assure to live a healthy life with style.

YOUR TOTAL HEALTHSTYLE PROVIDER



85
ACCURATE PRODUCTS



90
COUNTRIES WORLDWIDE



30
UNIQUE TECHNOLOGIES

www.rossmax.com/home



BLOOD PRESSURE MONITOR




PARR detects pulse arrhythmia, including Atrial Fibrillation (AFib), Premature Contractions (PC), Tachycardia (TACH), and Bradycardia (BRAD).

Detects

- AFib** Atrial Fibrillation
- PC** Premature Contraction
- TACH** Tachycardia
- BRAD** Bradycardia
- ARR** Non type - specific arrhythmia



THERMOMETER | PULSE OXIMETER



The continuous ear scanning mode views the entire ear canal and looks for the highest temperature that can be detected.



The non-contact temple forehead thermometer provides a consistent and longer distance temperature measurement.



ACT (Artery Check Technology) processes SpO₂ signals and determines the stiffness of the artery (Artery Condition) based on the derived wave form.

BLOOD GLUCOSE METER | NEBULIZER



The GDH-FAD test strips can eliminate interference from oxygen variations and non-glucose sugar such as maltose and galactose.



The VAT technology allows users to adjust different levels of nebulization rate at consistent particle size.



The 3BREATH is a unique 3in1 compressor nebulizer that integrates with a nasal washer and a nasal aspirator.

TECHNOLOGIES



Since 1988, Rossmax is committed to strive for the best technologies in each business field.

Providing value for the customer's specific need, is the key driver in inventing patented solutions. Rossmax does not compromise in fulfilling the highest quality standards and proofing the results through outstanding clinical validations.

www.rossmax.com/technologies



5IN1 TOTAL HEALTHCARE



BLOOD PRESSURE
MONITORING

SPO2 & ARTERY
MONITORING

WEIGHT & BODY FAT
MANAGEMENT

BLOOD GLUCOSE
MONITORING

TEMPERATURE
MEASUREMENT



Healthstyle APP



Health Station | 5in1

PREVENTION

The selective range of **FITNESS** products prevent chronic diseases such as **OBESITY** in the early stage.



Scale

MONITORING

The comprehensive product range in the field of **FEVER MANAGEMENT** and **HYPERTENSION** monitor health conditions on a professional and most accurate level.



Pulse Oximeter



Blood Pressure Monitor



Blood Glucose Meter



Thermometer



Calibration

PORTFOLIO

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Rossmax International Ltd. is under license. Other trademarks and trade names are those of their respective owners.

www.rossmax.com/products



THERAPY

Most effective solutions for **RESPIRATORY CARE** and **PAIN RELIEF**, complete the Rossmax Therapy range.



Nebulizer



Power Suction Unit



3in1 Respiratory Solution

PROFESSIONAL

A lean but sophisticated range of professional products including new initiatives for the **"LONG LIFE CARE"** category.



Blood Pressure Monitor



Aero Spacer




Peak Flow Meter



Accessories



Heating Pad



Air Mattress



Sphygmomanometer

KIDSHEALTH

A complete product range covering baby scale, nebulizer, ear and standard pen-type thermometers provide a total solution to **manage** the health of kids and babies easily.



Thermometer

Scale



Nebulizer

PORTFOLIO

www.rossmax.com/products





Rossmax' in-house Marketing agency is offering customized services to individual needs.
A very experienced, multi-lingual team of graphical designers, communication specialists, webmasters, movie-makers and social media experts are offering fast, creative and professional support.

design@rossmax.com







AFTER-SALES

In Rossmax, a sale is just the first step, but not the last. By providing excellent after-sales service we are committed to build a long-term relationship. During the last decades Rossmax has developed a dynamic global after-sales concept which consist on selling techniques, follow-up contact and complaint handling. Through education and regular feedbacks, the quality in each process has reached a very high level.




For more details, please contact your sales manager



VALIDATION

European Society of Hypertension & British Hypertension Society

Internal and local validations as well as proofs from well recognized institutions (such as BHS, ESH) confirm the accuracy, safety, stability, reliability and high performance of the entire product range. The Rossmax blood pressure monitors can be fully recommended for an application in the pregnant population too.



www.bihsoc.org

www.dablededucational.org

CONTACT US

SWITZERLAND

Rossmax Swiss GmbH
Widnauerstrasse 1
CH-9435 Heerbrugg, Switzerland

Tel: +41 71 544 43 25
emea@rossmax.com

USA

Rossmax USA
11545 Sorrento Valley Rd, Suite 308,
San Diego, CA. 92121

america@rossmax.com

TAIWAN

Rossmax International Ltd
12F, No. 189, Kang Chien Road
Taipei, 114, Taiwan

Tel: +886 2 2659 7888
Fax: +886 2 2659 7666
info@rossmax.com

JAPAN

Rossmax Japan
962-7 Shimoyasumatsu, Tokorozawa,
Saitama 359-0024 Japan

Tel: +81 042 9444 756
japan@rossmax.com

INDIA

Rossmax India Sub Continent
74, Lord Buddha Society-2nd Floor, Sector 21-C
Opp-Gymkhana Club-Faridabad 121003
NCR of Delhi, Haryana, India

Tel: +91 98116 00264
india@rossmax.com

CHINA

GMC Inc.
No. 686, Su Chu Rd., Chuzhou,
Anhui, China

Tel: +86 550 2172888
Fax: +86 550 2172889
info@rossmax.com

Want to discover more? Just visit us @



