

Super-cosy, high-temperature Heating Pad

Model: HP3040A (regular size)



- Overheating protection
- Quick heat-up
- 4 temperature settings
- Even heat distribution
- Cosy, removable and washable cover
- Automatic switch-off after 90 mins
- Ergonomic controller with temperature level display and LED indicator
- Regular size 30 x 40 cm
- CE approval



Overheating protection



Quick heat-up



4 temperature settings



Even heat distribution



Cosy, removable and washable cover



Automatic switch-off after 90 mins

Rossmax Group

rossmax
TAIWAN

rossmax
SWISS

rossmax
INDIA

rossmax
JAPAN

rossmax
AMERICA

rossmax 脉博士
CHINA

Rossmax International Ltd.
12F, No.189, Kang Chien Rd.,
Taipei 114, Taiwan.
Tel +886 2 2659 7888
Fax +886 2 2659 7666
info@rossmax.com

Rossmax Swiss GmbH
Tramstrasse 16, 9442 Berneck
Switzerland
Tel +41 71 747 11 93
Fax +41 71 747 11 91
emea@rossmax.com

www.rossmax.com

OBM_OP_Therapy_EN_ver1512




Heating Pad

Heat Therapy

Heat has been used as an effective way of therapy for long time. While the overall quality of warmth is related to comfort and relaxation, heat therapy provides both pain relief and healing benefits for many types of pain. Heat therapy can be provided quickly and inexpensively in the form of heating pads, which can be applied safely to the affected areas.


Comfortable, safe and effective

Rossmax Heating Pads offer constant and selective heat levels. Thanks to the "Quick heat-up" feature and the "Three Safety Levels", Rossmax Heating Pads deliver various therapy options. A Heating Pad that really gets warm easily and comfortably.



Three Safety Levels

- 1 Short circuit safety
- 2 Overheat protection via thermal fuse
- 3 PTC temperature monitoring



Heating Setting


4 heating levels can be selected and are displayed in the digital LCD.

Auto-off LED Indicator

A LED indicates the "Auto-Off" function is ON.

Cosy Fabric Cover

The cover offers luxurious softness and soothing comfort and is easy washable.



Super-cosy, high-temperature Heating Pad

Model: HP4060A (XXL size)



- Overheating protection
- Quick heat-up
- 4 temperature settings
- Even heat distribution
- Cosy, removable and washable cover
- Automatic switch-off after 90 mins
- Ergonomic controller with temperature level display and LED indicator
- XXL size 40 x 60 cm
- CE approval



Overheating protection



Quick heat-up



4 temperature settings



Even heat distribution



Cosy, removable and washable cover



Automatic switch-off after 90 mins