



HEALTH STATION	STL 501F (Deluxe)	STM 300C (Desktop)	STM 300F (Desktop)
	Hardware Specifications		
Size	(L) 412 x (D) 201 x (H) 864 mm	(L) 430 x (D) 310.5 x (H) 500 mm	
OS	Android 8.1		
CPU	1.5GHz		
RAM	2GB		
HHD	eMMC 16GB		
Network	LAN (Intel) x 1 - WLAN RTL8723BS 802.11 b/g/n 1T1R		
Bluetooth	4.0		
Speaker	2W		
Touch Screen	27" FHD Panel	15.6" FHD Panel	15.6" FHD Panel
Printer	●	●	●
4in1RACK	●	●	●
Camera	●	●	●
NFC Reader	●		●
Smart Card Reader	●		●
QR Code / Bar Code Scanner	●		●
	Health Devices		
X5 BT - Blood Pressure Monitor	●	●	●
HS200 - Blood Glucose Meter	●	●	●
SB210 - Pulse Oximeter	●	●	●
HC700 - Thermometer	●	●	●
WF262 - Body Fat Monitor	●	●	●
	Software Specifications - Stand alone		
Slideshow	●	●	●
Health Check Up	Blood Pressure, Blood Glucose, SpO ₂ , Temperature, Weight, Body & Visceral Fat, BMI, BMR, Body Muscle		
Print Receipt	Check up result, Coupon		
Settings	Load slides, Set coupon content, Change banner		



MEDICAL KIT	STP 100C
	Licensed Android System Software included
Recommended Tablet	Samsung Galaxy Tab A 8" (not included)
	Recommended Hardware Specifications for Tablets
OS	Android 9.0 above
CPU	2GHz
RAM	2GB
HHD	32GB
Network	WiFi 802.11 b/g/n
Bluetooth	4.2
Touch Screen	8" Panel (800*1280)
Camera	For reading Barcodes
	Health Devices included
X5 BT - Blood Pressure Monitor	●
HS200 - Blood Glucose Meter	●
SB210 - Pulse Oximeter	●
HC700 - Thermometer	●
WF262 - Body Fat Monitor	●
	Software Specifications - Stand alone
Health Check Up	Blood Pressure, Blood Glucose, SpO ₂ , Temperature, Weight, Body & Visceral Fat, BMI, BMR, Body Muscle
Check Up Results	By scanning the QR Code on the Tablet screen