

Car Charger User's Guide


Technical Specification

Input	DC 12V 1A
Output	DC 12V 1A
Operating Temperature	50°F to 104°F (10°C to 40°C)
Storage Temperature	-4°F to 140°F (-20°C to 60°C)
Size	94.3 x 19.7 x 22.6 mm
Weight	60g

Usage tips

Please refer to the following condition and the steps as usage tips.

1. To protect your device against sudden surges in current, power on the device only after the car engine has been started.
2. It's suggested to run the charging process while the car engine is turned on.



Cautionary Note

1. Only use the car charger to the authorized nebulizer models. Any unauthorized device is prohibited to connect with this car charger.
2. When switch off the car engine, please unplug the car charger from cigarette lighter jack.
3. Do not expose car charger to rain or excessive moisture.
4. Do not modify the car charger.
5. Do not cover the car charger which will cause risk of overheating.
6. Dispose of device according to local regulations.

Manufacturer: Dongguan Liwei Electronics Co., Ltd.
Room 5013, Floor 5th, Building . 1st Xihu Zhong Road,
Shilong Town, Dongguan City, Guangdong Province, P.R. C.



CarCharger@icubike_w@163